



DAI[®]

The Gold Standard

Volume 

Viscosity 

Device Type

Disposable Auto Injector

Primary Container

1mL Long Syringe with Staked Needle

Dosage Type

Fixed | Full Contents

Route of Administration

Subcutaneous

Needle Insertion

Automatic

Activation Method

Button Activated

Injection Feedback

Audible click at beginning & end of injection, viewing window

Handling Instructions



1
Remove cap



2
Hold device against injection site



3
Push & release activation button
to activate



+
Discard in sharps container

Key Features

- ✓ Simple 3-step operation
- ✓ Button activation
- ✓ Automatic needle insertion
- ✓ Audible, visual & tactile feedback
- ✓ One-handed operation
- ✓ Prevents needle-related injuries
- ✓ Permanently hidden needle
- ✓ Passive locking needle cover