



Molly®
2.25 mL
The ideal safety solution

Volume 

Viscosity 

Device Type

Disposable Auto Injector

Primary Container

2.25mL Long Pre-Filled Syringe | RNS Compatible

Dosage Type

Fixed | Full Contents

Route of Administration

Subcutaneous

Needle Insertion

Manual

Activation Method

Needle Cover

Injection Feedback

Audible click at beginning & end of injection, viewing window

Handling Instructions



Remove cap



Press & hold against injection site to activate



Discard in sharps container

Key Features

- ✓ Simple 2-step operation
- ✓ Needle cover activation
- ✓ Manual needle insertion
- ✓ Audible, visual & tactile feedback
- ✓ One-handed operation
- ✓ Larger volume
- ✓ Permanently hidden needle
- ✓ Ergonomic pull- or twist-off cap
- ✓ Robust anti-roll design