

**DAI®**  
*The Gold Standard*

 **SHL GROUP**  
www.shl-group.com



## *Technical Specifications*

|                          |   |
|--------------------------|---|
| <b>Device Type</b>       | Disposable Autoinjector                       |
| <b>Primary Container</b> | 1 mL Long Pre-filled Syringe   RNS Compatible |
| <b>Dosage Type</b>       | Fixed   Full Contents                         |
| <b>Needle Extension</b>  | Up to 1/2" or 12.7mm   Subcutaneous           |
| <b>Needle Insertion</b>  | Automatic                                     |
| <b>Needle Removal</b>    | Manual  |
| <b>Activation Method</b> | Button  |

Copyright © 2015 SHL Group. All rights reserved. Version 1.1 Information accurate at time of print and is subject to change without notice. Available as standard or customized to client specifications. Contact [info@shl-group.com](mailto:info@shl-group.com) for more information.

# Handling Instructions



Remove cap



Press & hold against injection site



Push & release activation button while holding device on injection site



Remove device from injection site & discard in sharps container

## Key Features

- ✓ Permanently hidden needle
- ✓ Portable design
- ✓ Audible, visual & tactile feedback
- ✓ Compatible with Rigid Needle Shields