

# Molly®

*The Ideal Safety Solution*



## *Technical Specifications*

<b>Device Type</b>	Disposable Autoinjector
<b>Primary Container</b>	1mL Long Pre-filled Syringe   RNS Compatible
<b>Dosage Type</b>	Fixed   Full Contents
<b>Needle Extension</b>	Up to 1/4" or 6mm   Subcutaneous
<b>Needle Insertion</b>	Manual
<b>Needle Removal</b>	Manual
<b>Activation Method</b>	Needle cover

# Handling Instructions



1.

Remove cap



2.

To activate, press  
needle cover against  
injection site & hold



+

Remove device  
from injection site  
& discard in sharps  
container

## Key Features

- ✓ Ultra compact design
- ✓ Needle cover activation
- ✓ Pre-configured or customized to client specifications
- ✓ Compatible with Rigid Needle Shields