

SDI MIX[®]+NIT[™]

Twist & Mix



Technical Specifications

Device Type	Disposable Autoinjector
Primary Container	1mL Dual Chamber Cartridge
Dosage Type	Fixed Full Contents
Needle Extension	Up to 1/2" or 12.7mm Subcutaneous
Needle Insertion	Manual
Needle Removal	Manual
Activation Method	Needle Cover

Handling Instructions



1. Turn rear dial to mix



2. Twist needle cap off.
Needle cover will
automatically extend.
Device will prime.



3. Press & hold against
injection site



+
Remove device
from injection site
& discard in sharps
container

Key Features

- ✓ Needle Isolation Technology [NIT™]
- ✓ Automatic reconstitution & priming
- ✓ Needle cover activation
- ✓ Audible, visual & tactile feedback